

SPECIFICATIONS

- ▶ Chemical formula: $ZrO_2 / Y_2O_3 / Al_2O_3$
- ▶ Chemical name: Yttrium stabilized tetragonal zirconium dioxide
- ▶ Description of product: Presintered blanks made of zirconium dioxide for dental use Typ II, Class 4 acc. to EN ISO 6872

CHEMICAL CHARACTERISTICS (TYP.)

$ZrO_2 / HfO_2 / Y_2O_3$	[wt%]	> 99.0
Y_2O_3	[wt%]	9.30 ± 0.30
HfO_2	[wt%]	< 5.0
Al_2O_3	[wt%]	< 0.1
Fe_2O_3	[wt%]	< 0.01

MECHANICAL & PHYSICAL CHARACTERISTICS (TYP.)

Density	[g/cm ³]	> 6.0
4-point bending strength	[MPa]	600 ± 100
Young's modulus	[Gpa]	210
Thermal expansion	[10 ⁻⁶ 1/K]	10.4
Open porosity		none
Radio activity		
Activity ²³⁸ U	[Bq/g]	< 1.0
Biocompatibility and chemical solubility		
Cytotoxicity ISO ¹		Inoffensive
Chemical analysis ²	[µg/cm ²]	< 2,000

¹ (DIN EN ISO 10993 / DIN EN ISO 9363)

² (DIN EN ISO 6872)

PICTURE DENTAL BLANKS



PRODUCT DESCRIPTION & APPLICATIONS

Presintered blanks from tetragonal stabilized zirconium oxide are a semi-finished ceramic product with special properties for medical applications, especially for machining of crowns and bridges with CAD/CAM milling systems. They are suitable for the fabrication of crowns and bridges up to 3 units in monolithic version or as framework for coating with ceramic. The high translucency and the good biocompatibility are outstanding characteristics and support the suitability for medical applications.



STARCERAM® Z-SMILE

2023.01.24

PACKAGING

StarCeram® Z-Smile products are individually packed, depending on size.
Packaging acc. to transport regulations.

STORAGE & HANDLING

Handle with care! Ceramic products may break if not treated properly. Storage and handling are subject to the rules and regulations in the country of use.

HAZARDS IDENTIFICATION IN ADVERTISING (REGULATION (EC) NO 1272/2008 ARTICLE 48)

none

DOCUMENTATION

An inspection document in accordance with EN 10204 is supplied with every shipment.

The information contained herein, and in particular the recommendations relating to the application and end-use of our products, are given in good faith based on our current knowledge and experience. The information do not constitute a guarantee with respect to properties of the products. Since we have no control over the application modalities, no guarantee is granted with respect to merchantability or fitness for a particular purpose. It is your sole responsibility to validate the suitability and completeness for your own use. We therefore recommend to always perform a test according to specific circumstances. Any liability in respect of the information in this data sheet or any other written or oral recommendation(s) regarding the concerned product is excluded, except if otherwise explicitly agreed and except in case of death or personal injury, intent or gross negligence and under any applicable mandatory product liability law. Technical data are subject to change without notice. The newest version of the technical data sheet replaces all preceding versions. The trademarks, trade names, logos and other designations of origin contained in this data sheet are registered and unregistered intellectual property rights of KYOCERA Fineceramics Europe GmbH. It is forbidden to copy or use information from this data sheet in whole or in part, especially in dealings with third parties.

KYOCERA Fineceramics Europe GmbH · Standort/Plant Selb · 95100 Selb / Germany
Lorenz-Hutschenreuther-Strasse 81 · Tel.: +49 (0)9287 - 807-0 · E-Mail: info@kyocera-precision.com
www.kyocera-precision.com

© 2023 KYOCERA Fineceramics Europe GmbH. We reserve the right to make changes without prior notice.